

## IMMUNOFLUORESCENCE PROTOCOL

### A. Cultured Cells

1. Coat coverslips with polyethyleneimine or poly-L-lysine for 1 hr at room temperature.
2. Rinse coverslips well with sterile H<sub>2</sub>O (3 times 5 min each).
3. Allow coverslips to dry completely and sterilize them under UV light for at least 4 hrs.
4. Grow cells on glass coverslips or prepare cyospin or smear preparation.
5. Rinse briefly in phosphate-buffered saline (PBS).

### B. Fixation

1. Fix the samples either in ice-cold methanol, acetone (1-10 min) or in 3-4% paraformaldehyde in PBS pH 7.4 for 15 min at room temperature.
2. Wash the samples twice with ice cold PBS.

### C. Permeabilization

If the target protein is expressed intracellularly, it is very important to permeabilize the cells. Note: acetone fixed samples do not require permeabilization.

1. Incubate the samples for 10 min with PBS containing 0.25% Triton X-100 (or 100  $\mu$ M digitonin or 0.5% saponin). Triton X-100 is the most popular detergent for improving the penetration of the antibody. However, it is not appropriate for the use of membrane-associated antigens since it destroys membranes.
2. Wash cells in PBS three times for 5 min.

### D. Blocking and Incubation

1. Incubate cells with 1% BSA in PBST for 30 min to block unspecific binding of the antibodies (alternative blocking solutions are 1% gelatin or 10% serum from the species that the secondary antibody was raised in).
2. Incubate cells in the diluted antibody in 1% BSA in PBST in a humidified chamber for 1 hour at room temperature or overnight at 4°C.
3. Decant the solution and wash the cells three times in PBS, 5 min each wash.
4. Incubate cells with the fluorochrome-conjugate secondary antibody in 1% BSA for 1 hour at room

temperature in dark.

5. Decant the fluorochrome-conjugate secondary antibody solution and wash three times with PBS for 5 min each in dark.

#### **E. Counter staining**

1. Incubate cells on 0.1-1  $\mu\text{g/ml}$  Hoechst or DAPI (DNA stain) for 1 min.
2. Rinse with PBS.

#### **F. Mounting**

1. Mount coverslip with a drop of mounting medium.
2. Seal coverslip with nail polish to prevent drying and movement under microscope.
3. Store in dark at  $-20$  or  $4^{\circ}\text{C}$ .
4. Test under fluorescence microscope.