

## FLOW CYTOMETRY PROTOCOL

Indirect flow cytometry requires two incubation steps, firstly with a primary antibody then with a compatible secondary antibody. The secondary antibody has the fluorescent dye (FITC, PE, Cy5, etc.) conjugated.

### A. Fixation

1. Harvest and wash the cells then determine the total cell number.
2. Resuspend the cells to approximately  $1-5 \times 10^6$  cells/ml in ice cold PBS.
3. Resuspend cells briefly in 0.5-1 ml PBS. Add formaldehyde to a final concentration of 2-4% formaldehyde.
4. Fix for 10 minutes at 37°C.
5. Chill tubes on ice for 1 minute.

**Note:** For extracellular staining with antibodies that do not require permeabilization; for intracellular staining, proceed to permeabilization.

### B. Permeabilization

1. Permeabilize cells by adding ice-cold 100% methanol slowly to pre-chilled cells, while gently vortexing, to a final concentration of 90% methanol. Alternatively, to remove fix prior to permeabilization, pellet cells by centrifugation and resuspend in 90% methanol.
2. Incubate 30 minutes on ice.
3. Proceed with staining or store cells at -20°C in 90% methanol.

### C. Immunostaining

1. Add 0.1-10 µg/ml of the primary antibody. Dilutions, if necessary, should be made in 3% BSA/PBS.
2. Incubate for at least 30 min at room temperature or 4°C in the dark.
3. Wash the cells 3-times by centrifugation at 400 g for 5 min and resuspend them in ice cold PBS. You may need to adjust the conditions of the centrifugation for the cell types used.
4. Dilute the fluorochrome-labeled secondary antibody in 3% BSA/PBS at the optimal dilution and then resuspend the cells in this solution.
5. Incubate for at least 20-30 minutes at room temperature of 4°C. This incubation must be done in the dark.
6. Wash the cells 3 X by centrifugation at 400 g for 5 min and resuspend them in ice cold PBS, 3% BSA, 1% sodium azide.
7. Store the cell suspension immediately at 4°C in the dark.



WEBSITE  
E-MAIL

[www.immunoway.com](http://www.immunoway.com)  
[order@immunoway.com](mailto:order@immunoway.com)  
[tech@immunoway.com](mailto:tech@immunoway.com)

#### **D. Optional DNA Stain**

1. Resuspend cells in 0.5 ml of DNA dye.
2. Incubate for at least 5 minutes at room temperature.
3. Analyze cells in DNA stain on flow cytometer.