**Applications** 

IF,FC,WB



# Rabbit Anti Dog IgG(H+L) (AbFluor 680)

CatalogNo: RS4015

#### **Key Features**

Host Species Reactivity

RabbitDog

Isotype Conjugate
• IgG • AbFluor 680

#### **Recommended Dilution Ratios**

Suggested starting 1:50-1:1000 dilutions for most fluorescent applications while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.

### Storage

**Storage\*** Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of

product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot

to avoid repeated freezing and thawing

**Formulation** 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50%

glycerol, 0.02% Sodium Azide

#### **I** Basic Information


**Clonality** Polyclonal

## Immunogen Information

### | Target Information

**Protein Name** 

#### **Validation Data**



To use the Alexa Fluors with fluorescent imagers, use a spectral line of the blue laser diode for Alexa Fluors 405, a cyan (488 nm) laser for Alexa Fluors 488, a yellow (526 nm) laser for Alexa Fluor 550 or 594, and a red (633 nm) laser for Alexa Fluor 649. The Alexa Fluor 680 and 790 fluors are compatible with laser- and filter-based infrared imaging instruments that emit in the 700 nm, and 800 nm

#### **Contact information**

Orders: order@immunoway.com
Support: tech@immunoway.com

Telephone: 877-594-3616 (Toll Free), 408-747-0185

Website: http://www.immunoway.com

Address: 2200 Ringwood Ave San Jose, CA 95131 USA



Please scan the QR code to access additional product information:

Rabbit Anti Dog
IgG(H+L) (AbFluor
680)

For Research Use Only. Not for Use in Diagnostic Procedures.

Antibody | ELISA Kits | Protein | Reagents