Goat Anti Rabbit IgG(H+L) (AbFluor 594)

CatalogNo: RS3611

Key Features

Host Species
- Goat

Reactivity
- Rabbit

Applications
- IF,ELISA,FC

Conjugate
- AbFluor 594

Storage

Storage*	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Recommended Dilution Ratios

IF 1:200-1:1000
Flow Cyt 1:100-1:1000
ELISA (Use at an assay dependent concentration)


Basic Information

Immunogen Information

Target Information

Protein Name

Validation Data



Inhaled Lipid Nanoparticles Alleviate Established Pulmonary Fibrosis Small Yang Yang IF Mouse murine lung epithelial-12 (MLE-12) cell

AbFluor	350	346/442	Blue
AbFluor	405	401/421	Blue
AbFluor	488	496/519	Green
AbFluor	532	532/553	Yellow
AbFluor	555	555/565	Yellow
AbFluor	568	578/603	Red/Orange
AbFluor	594	590/617	Red/Orange
AbFluor	633	632/647	Red
AbFluor	647	650/665	Red
AbFluor	660	663/690	Near IR
AbFluor	680	679/702	Near IR
AbFluor	750	749/775	Near IR
AbFluor	790	784/814	Near IR

To use the AbFluors with fluorescent imagers, use a spectral line of the blue laser diode for AbFluors 405, a cyan (488 nm) laser for AbFluors 488, a yellow (526 nm) laser for AbFluor 550 or 594, and a red (633 nm) laser for AbFluor 649. The AbFluor 680 and 790 fluors are compatible with laser- and filter-based infrared imaging instruments that emit in the 700 nm, and 800 nm

Contact information

Orders: order@immunoway.com
Support: tech@immunoway.com
Telephone: 877-594-3616 (Toll Free), 408-747-0185
Website: http://www.immunoway.com
Address: 2200 Ringwood Ave San Jose, CA 95131 USA

Please scan the QR code to access additional product information: **Goat Anti Rabbit IgG(H+L) (AbFluor 594)**

For Research Use Only. Not for Use in Diagnostic Procedures.